

# Great BIONUTRON contest, Win a 1-year supply of BIONUTRON !



The winner will be drawn from among all correct answers received by 30 April 2010.  
The winner may choose any one of the 6 existing Bionutron barrel references (of a maximal 400-euro value).

Send your answers:

- by email at [info@audevard.com](mailto:info@audevard.com) (mentioning the question number and the chosen answer).
- by fax at : +33 1 47 56 38 39.
- by post : Audevard labs, 42 rue Mederic, 92 110 Clichy, FRANCE.

Yes, I would like to receive information from AUDEVARD Labs (news, products, offers, etc.)

Last name \_\_\_\_\_

First name \_\_\_\_\_

Email \_\_\_\_\_

Country \_\_\_\_\_

## 1 Why use a supplement from the BIONUTRON range?

- A.  To balance the ration and optimize the horse's well-being
- B.  To replace the food
- C.  To excite the horse

## 2 How are BIONUTRON products manufactured?

- A.  By heating raw materials
- B.  By means of a pasteurisation process
- C.  By means of a cold process to preserve the vitamins

## 3 What is the purpose of the nutrigenomic ingredients contained in BIONUTRON supplements?

- A.  To increase the fertility of the mare or stallion
- B.  To optimize the utilization of the horse's genetic potential
- C.  To draw up the horse's genetic map



This information will be used exclusively by Audevard Labs.  
For further information, you can contact AUDEVARD: [info@audevard.com](mailto:info@audevard.com).  
Winners will be drawn from all correct answers and notified by email.